

---

## CERTIFICATE OF ANALYSIS


**Product: DMEM-LOW GLICOSE – 1 X Concentrated  
With Sodium Bicarbonate (3.7 g/L)**

Lot # 040.717TJA

Catalog # BR30002.05

Manufacture Date: Jun/29/2017  
Expiration Date: Jun/29/2018

Test	Specification	Units	Results
Appearance	Clear, transparent		Clear, transparent
pH w/ NaHCO <sub>3</sub>	7.2 +/- 0.2		7.20
Bacteria and yeast	Negative		Negative
Endotoxins	≤ 10.0	EU/mL	≤ 2,0
<i>Mycoplasma</i>	Negative		Negative
Osmolality w/NaHCO <sub>3</sub>	305 +/- 5%	mOsm/Kg H <sub>2</sub> O	310



Responsável Técnico  
LGC Biotecnologia Ltda.