



---

## CERTIFICATE OF ANALYSIS

**Product: Phosphate Saline Buffer Dulbecco's Modification,  
Without Calcium and Magnesium**

**Lot # 290.417BJ**

**Catalog # BR30013-05**

**Manufacture Date: 2017/Abr/24**

**Expiration Date: 2018/Abr/24**

<b>Test</b>	<b>Specification</b>	<b>Units</b>	<b>Results</b>
Appearance	Clear, transparent		Clear, transparent
pH w/o NaHCO <sub>3</sub>	7.2 +/- 0.2		7.21
Bacteria and yeast	Negative		Negative
Endotoxins	≤ 10.0	EU/mL	Negative
<i>Mycoplasma</i>	Negative		Negative
Osmolality	270 +/- 10%	mOsm/Kg H <sub>2</sub> O	280

---

Assinado eletronicamente